

William A. Rutala, PhD, MPH,^{1,2} **Hajime Kanamori**, MD, PhD, MPH,^{1,2} Maria F. Gergen, MT (ASCP),¹ Emily E. Sickbert-Bennett, PhD,^{1,2}

Daniel J. Sexton, MD,³ Deverick J. Anderson, MD, MPH,³ David J. Weber, MD, MPH,^{1,2} and the CDC Prevention Epicenters Program

¹Hospital Epidemiology, University of North Carolina Health Care, Chapel Hill, NC, ²Division of Infectious Diseases, University of North Carolina School of Medicine, Chapel Hill, NC

³Duke Infection Control Outreach Network, Division of Infectious Diseases, Duke University Medical Center, Durham, NC, USA

William A. Rutala, Ph.D., M.P.H.
Director, Hospital Epidemiology,
Occupational Health, and Safety Program
Room 1001 West Wing, Manning Drive, UNC Hospitals
Chapel Hill, NC 27514
Phone: 919-843-1397 Email: brutala@unch.unc.edu



THE UNIVERSITY
of NORTH CAROLINA
at CHAPEL HILL

Background

- Recent evidence has demonstrated that the contaminated hospital environment is an important source for transmission of epidemiologically-important pathogens (EIP) such as MRSA and VRE.
- An overhead light fixture technology, which continuously and safely disinfects the environment using light-emitting diodes (LEDs) that emit a high-intensity, narrow spectrum (HINS) light has been proposed as an infection prevention strategy.
- This technology uses LEDs to create a narrow bandwidth of high-intensity visible violet light with a peak output of 405nm that react with porphyrin molecules to generate reactive oxygen species that kill microorganisms.
- The purpose of this evaluation was to determine the effectiveness of HINS-light for the reduction of epidemiologically-important pathogens.

Methods

- The new technology was evaluated in two different clinical configurations (White Disinfection Mode and Supplemental Blue-Disinfection Mode).
- In Phase 1, two 2x2 blended-white, ceiling-mounted fixtures were used which provided both disinfection and ambient white illumination for use in normal clinical conditions in an occupied room (surface irradiance ~0.12-0.16 mW/cm² measured at the pathogen location).
- In Phase 2, a higher-level of disinfection was studied with the addition of a 2x4 overhead. This fixture emits only disinfecting (blue) light (surface irradiance=0.34-0.44 mW/cm²).
- The four test organisms were *C. difficile* spores (BI), MRSA (ATCC 43300), VRE (ATCC 51299), and MDR-*A. baumannii* (MDRA). Formica test surfaces were inoculated with 100-500 CFUs of test organisms.
- Once dry, triplicate samples were collected with Rodac plates containing DE Neutralizing Agar at times 0, 1hr, 3hr, 5hr, 6hr, 7hr, 24hr, 48hr, and 72hr then appropriately incubated.
- To accommodate the natural die-off of vegetative bacteria, a control test surface (Formica) was placed in an adjacent area but not exposed to the light disinfection.

Results




Figure 1. Survival in MDR-*Acinetobacter* (A), MRSA (B), VRE (C), and *C. difficile* (D) with Blue Light

Table 1. Time to specified percent reductions of epidemiologically-important pathogens with “blue” light and “white” light.

Treatment	Pathogen	Time (least number of hours) to achieve microbial reduction				Maximum reduction achieved (%)
		25%	50%	90%	100%	
Blue light	MRSA		3	48	48	100
	VRE	5	24	24	48	100
	MDR- <i>Acinetobacter</i>	1	5	NA	NA	88
	<i>C. difficile</i>	5	72	NA	NA	65
White light	MRSA	7	24	48	72	100
	VRE	24	NA	NA	NA	47
	MDR- <i>Acinetobacter</i>	6	24	48	72	100
	<i>C. difficile</i>	5	NA	NA	NA	25

NA, not achieved. MRSA, methicillin-resistant *Staphylococcus aureus*; VRE, vancomycin-resistant *Enterococcus*, MDR-*Acinetobacter*, multidrug resistant *Acinetobacter*.




Figure 2. Survival in MDR-*Acinetobacter* (A), MRSA (B), VRE (C), and *C. difficile* (D) with White Light

Results Summary

- These results demonstrated that the 405nm light inactivated three vegetative bacteria (MRSA, VRE, MDRA) on surfaces with contact times of 1-72hr.
- Statistical differences ($p < 0.05$) were observed using “blue” light for VRE at 24 and 48 hours, for MRSA at 3, 5, 6 and 7 hours, for MDR-*Acinetobacter* at 5, 6, 7, and 24 hours, and for *C. difficile* spores at 5 and 72 hours.
- The “white” light demonstrated no statistically improved reduction (i.e., $p > 0.05$) at any time point compared to control die-off for VRE, MDR-*Acinetobacter*, and *C. difficile* but there was a statistical difference for MRSA at 1, 7, 24 and 48 hours.
- The inactivation was more significant when the surface irradiance was increased by adding the blue light.

Conclusions

- We demonstrated that the “blue” light significantly reduced both vegetative bacteria and spores at some time points over a 72 hour exposure period.
- In addition to episodic disinfection (e.g., UV), this continuous light disinfection technology could be considered for several healthcare decontamination applications (e.g., OR).
- Given that environmental surfaces in patient’s room are often not thoroughly disinfected plus that recontamination occurs rapidly it is important to develop either methods of continuous disinfection or a germicide with persistent antimicrobial effectiveness.
- Whether these reductions are sufficient to reduce healthcare-associated infections requires further studies.

Acknowledgments: The authors thank Dr. Jeff Laux for his assistance with statistical analysis.

October 26-30 • New Orleans, LA • www.idweek.org