

Addressing pulmonary nocardiosis risk in immunocompromised patients: Development and validation of a commercially available PCR

Emily M. Smith, Katelyn Bartlett, James Grantham, Natalie Powell, Michelle Altrich, Steve Kleiboeker, Mark C. Wissel
Viracor Eurofins, Lee's Summit, MO USA



Eurofins
Clinical Diagnostics

INTRODUCTION

Current state of diagnosis of pulmonary nocardiosis. Pulmonary nocardiosis is an infection targeting immunocompromised patients that is characterized by high mortality and requires non-frontline antibiotics for treatment, such as linezolid and trimethoprim/sulfamethoxazole (Martinez T et al., 2007). Early diagnosis of *Nocardia* infections remains a major challenge. Difficulties persist in identifying *Nocardia* due to extended culture times for BAL.

Importance of early diagnosis for nocardiosis. Nocardiosis is currently confirmed or excluded by BAL fluid culture followed by further phenotypic identification steps. A culture-independent method with more timely results would accelerate the administration of appropriate treatment (Wilson JW, 2012). A rapid *Nocardia* PCR assay for BAL has neither been validated nor offered for clinical testing to our knowledge.

Early detection of nocardiosis aided by rapid molecular testing. This commercially available real-time PCR assay may permit early detection of nocardiosis and administration of appropriate treatment in a timely manner, when compared to the culture method.

MATERIALS AND METHODS

A rapid and sensitive real-time PCR assay was developed by Viracor Eurofins Clinical Diagnostics for detection of *Nocardia* species by TaqMan technology for use in BAL specimens. Oligonucleotides for a *Nocardia* PCR were aligned to the 16S regions of common causative agents of nocardiosis and over 30 others (Brown-Elliott BA et al., 2012). The target species include: *N. abscessus*, *N. africana*, *N. asteroides*, *N. brasiliensis*, *N. cyriacigeorgica*, *N. farcinica*, *N. nova*, *N. otitidiscaviarum*, *N. pseudobrasiliensis*, *N. transvalensis*, and *N. wallacei*. The assay has one fluorogenic TaqMan probe that is FAM-labeled and thus detects, but does not distinguish *Nocardia* species (Table 1). Rapid automated nucleic acid extraction (<1 hour for 24 samples) followed by fast PCR (<1 hour) was validated according to relevant compliance standards at Viracor Eurofins (Figure 1). Spiked and un-spiked human BAL samples were used to assess analytical specificity, limit of detection, precision and accuracy using commercially available *Nocardia* (Table 2) and non-*Nocardia* strains (Table 3).

Table 1: *Nocardia* PCR Assay Oligonucleotides

Assay Oligonucleotides	Number
<i>Nocardia</i> Forward Primer	1
<i>Nocardia</i> Reverse Primer	1
<i>Nocardia</i> Probes (FAM)	4

0.5mL BAL → Nucleic Acid → Detected / Not Detected

Figure 1. *Nocardia* Real-Time PCR Assay Workflow

RESULTS

Table 2. Analytical Specificity: Inclusivity

Detected	Detected in-Silico
<ul style="list-style-type: none"> <i>Nocardia abscessus</i> <i>Nocardia africana</i> <i>Nocardia asteroides</i> <i>Nocardia brasiliensis</i> <i>Nocardia cyriacigeorgica</i> <i>Nocardia farcinica</i> 	<ul style="list-style-type: none"> <i>Nocardia nova</i> <i>Nocardia otitidiscaviarum</i>¹ <i>Nocardia pseudobrasiliensis</i> <i>Nocardia transvalensis</i> <i>Nocardia wallacei</i>¹
	<ul style="list-style-type: none"> <i>Nocardia puris</i> <i>Nocardia cerradoensis</i> <i>Nocardia vulneris</i> <i>Nocardia brevicatena</i> <i>Nocardia teneriffensis</i> <i>Nocardia araoensis</i> <i>Nocardia testacea</i> <i>Nocardia arthritis</i> <i>Nocardia beijingensis</i> <i>Nocardia mexicana</i> <i>Nocardia amamiensis</i>
	<ul style="list-style-type: none"> <i>Nocardia crassostraea</i> <i>Nocardia thailandica</i> <i>Nocardia higoensis</i> <i>Nocardia harenae</i> <i>Nocardia carnea</i> <i>Nocardia asiatica</i> <i>Nocardia vinacea</i> <i>Nocardia aobensis</i>

¹ *N. otitidiscaviarum* and *N. wallacei* were detected below the LOD of the assay.

Table 3. Analytical Specificity: Cross-Reactivity

Detected ¹	Not Detected
<ul style="list-style-type: none"> <i>Crossiella cryophila</i> <i>Rhodococcus equi</i> 	<ul style="list-style-type: none"> <i>Acinetobacter baumannii</i> <i>Actinomyces israelii</i> <i>Arcanobacterium haemolyticum</i> <i>Bordetella holmesii</i> <i>Bordetella parapertussis</i> <i>Bordetella pertussis</i> <i>Chlamydia pneumoniae</i> <i>Chlamydia trachomatis</i> <i>Corynebacterium diphtheriae</i> <i>Fusobacterium necrophorum</i> <i>Haemophilus influenzae</i> Type B <i>Klebsiella pneumoniae</i> <i>Legionella pneumophila</i> <i>Moraxella catarrhalis</i> <i>Mycobacterium bovis</i> <i>Mycobacterium microti</i> <i>Mycobacterium tuberculosis</i>
	<ul style="list-style-type: none"> <i>Mycoplasma pneumoniae</i> <i>Neisseria gonorrhoeae</i> <i>Neisseria meningitidis</i> <i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Streptococcus dysgalactiae</i> <i>Streptococcus pneumoniae</i> <i>Streptococcus pyogenes</i> <i>Saccharothrix longispora</i> <i>Nocardiosis dassonvillei</i> subsp. <i>dassonvillei</i> <i>Nocardiosis dassonvillei</i> subsp. <i>albirubida</i>

¹ Two non-*Nocardia* species were detected at low levels when amounts of the organism were tested at the highest concentrations permitted by available bacterial stocks. *Rhodococcus equi* and *Crossiella cryophila* have had >100 and zero reported cases in literature, respectively.

Table 4. Limit of Detection

<i>Nocardia</i> species	Specimen type	Copies per mL BAL
<i>Nocardia cyriacigeorgica</i> ¹	BAL	206
<i>Nocardia nova</i> ²	BAL	41
<i>Nocardia transvalensis</i> ²	BAL	26

¹ The limit of detection for *N. cyriacigeorgica* was determined by testing 10-20 replicates of human BAL spiked at 9 levels of *N. cyriacigeorgica*. A total of 220 replicates were extracted and assayed by *Nocardia* PCR.

² A partial limit of detection for *N. nova* and *N. transvalensis* was performed, as these, along with *N. cyriacigeorgica*, are common causative agents of nocardiosis. The results were analyzed by Probit to predict the limit of detection, which is defined here as the concentration at which 95% of replicates are detected.

RESULTS

Table 5. Qualitative Analytical Cutoff

<i>Nocardia</i> species	Specimen type	Copies per mL BAL
<i>Nocardia cyriacigeorgica</i> ¹	BAL	206
<i>Nocardia nova</i>	BAL	41
<i>Nocardia transvalensis</i>	BAL	26

Note. Qualitative analytical cutoff was determined from Probit analyses from Limit of Detection data (Table 4). Qualitative analytical cutoff is calculated as the C₉₅ upper 95% confidence interval for each of the three species tested.

Table 6. Precision: Intra- and Inter-Assay Reproducibility

Precision	BAL	
	High Concentration (200,000 copies/mL)	Low Concentration (2,000 copies/mL)
Intra % CV	<10%	<8%
Inter % CV	<21%	<26%

Note. Five replicates of low and high concentrations of *N. cyriacigeorgica* spiked in human BAL were analyzed over three experimental days by multiple operators/instruments. The % coefficient of variation (%CV) was determined using the mean and standard deviation values of the copies/mL.

Table 7. Analytical Accuracy

Analytical Accuracy	BAL	
	Positive	Negative
Detected	100% (33/33)	0% (0/11)
Not Detected	0% (0/33)	100% (11/11)

Note. Target species listed in the methods were spiked in human BAL in triplicate, and 11 un-spiked samples were analyzed for accuracy.

CONCLUSIONS

The specificity, inclusivity, sensitivity, precision and accuracy of a qualitative Real-time *Nocardia* PCR have been deployed as an aid in the diagnosis of pulmonary nocardiosis.

Nocardia PCR allows for a culture-independent method that can rapidly detect clinically relevant *Nocardia* species with an improved turnaround time, leading to prompt diagnosis and administration of appropriate treatment.

References

- Martinez T et al. Pulmonary nocardiosis: risk factors and outcomes. *Respirology*, 2007 May; 12(3): 394-400.
- Wilson JW. Nocardiosis: updates and clinical overview. *Mayo Clin Proc*. 2012 Apr;87(4):403-7.
- Brown-Elliott BA, Biehle J, Conville PS, Cohen S, Saubolle M, Sussland D, Wengenack N, Kriel K, Bridge L, McNulty S, Vasireddy R, Wallace RJ Jr. Sulfonamide resistance in isolates of *Nocardia* spp. from a US multicenter survey. *J Clin Microbiol*. 2012 Mar;50(3):670-2.